Supporting Harness For Cattle

0704804

NL	Ondersteuningsharnas voor koeien	
	Gebruiksaanwijzing	2
EN	Supporting harness for cattle	
	Manual.....	4
DE	Unterstützungsharnas für Kühe	
	Gebrauchsanweisung.....	6
FR	Harnachement de soutien pour bovins	
	Mode d'emploi.....	8
IT	Fascia per sollevamento mucche	
	Manuale di istruzioni.....	10

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

Voorbeeld 2:

Supporting harness for cattle

MS supporting harness for cattle is unique in that it offers excellent support of the hind quarters. The hoist supports the cow along its length as well as between direction as well as between the hind legs and around the stomach.

The harness is used to hoist cows that have fallen or that are temporarily paralysed and to keep them standing upright over a longer period of time. Cow lifting hoists are also used for this purpose but they should not be used for long periods of time as this damages the hips. The harness is used not only by cattle breeders, but also by veterinarians and fire-fighters. In the event of the above situation, fast and effective action is essential in order to achieve good and quick recovery. This is why no dairy farm should be without a cow harness.

The hoist supports the cow along its length as well as between the hind legs and around the stomach. The udder is lifted in the harness during hoisting, which spreads the weight evenly over the full length of the cow.

This even spread of weight means that the stomach does not have to carry the full weight during hoisting, which could otherwise put pressure on the heart and lungs. The cow is also lifted with a straight back and the udder is not allowed to hang. The special shape of the harness places very little pressure on the udder, even when the cow has recently calved.

During hoisting the harness supports in such a way that the straps do not cut in to the sides of the animal, thereby avoiding the blood flow from being cut off.

Milking or sucking of a calf

Milking or suckling of a calf is not a problem as the back flap, which covers the udder, can be untied and let down. The cow can be held upright for a long enough period of time while the back flap is down to allow the animal to be milked.

The back flap should be rolled to the front during milking or suckling. The long back straps can be buckled from the hind quarters of the cow to the rings, which stops the cow from sliding to the back.

The harness is made of heavy woven PVC fibres, that are reinforced with nylon straps. This material ensures a large load capacity so that cows of all breeds and weights can be easily hoisted.

The special construction of the harness ensures complete support of the cow during transport. The cow is comfortable and the pressure on heart and lungs is low.


Applying the harness

Make a firm hoisting construction on, for instance, the front loader of a tractor. The construction must be of sufficient strength for instance a metal tube beam. The length varies can be varied from 1.20 to 1.60 metres depending on the length of the cow. The harness can be put in the correct position under the cow by rolling her over and by pulling the harness underneath her. A cow lift can be used to hoist the cow for a short period of time.

1. The harness is fitted with 4 hoisting rings (2 rings at the height of the front legs and 2 rings at the height of the hind legs). These hoisting rings are to be attached to the hoisting construction at points A and B respectively by means of 2 ropes with hoisting hooks. The outside of the harness is reinforced with nylon straps.

2. Lead the 2 straps at the front between the front legs. Attach the hooks at point A of the hoisting construction. Use the spanner to adjust the strap to the correct length and tension.

3. Lead the 2 straps at the back between the hind legs. Cross the straps over each other to give an even distribution of weight. Make sure that the blood flow in the udder and the tail are not being cut off by the hoisting straps. Attach the hooks at point C of the hoisting construction. Use the spanner to adjust the strap to the correct length and tension.

Adjust all straps and hoisting ropes in such a way that the cow hangs horizontally. The cow can now be lifted upwards.

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**Schippers Italia Srl**
**Schippers Schweiz GmbH**
**Schippers UK Ltd.**
**LLC Schippers Russia**
**Schippers Canada Ltd**
7102 52nd Street Bay # 18 • Lacombe, AB • T4L 1Y9