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MS Gas operated tail docker standard

A complete docking instrument for quick and painless docking of newborn piglets

- Docks and cauterizes in a single operation, reducing the chance of infection.
- Operates on regular lighter fluid.
- Economical: docks up to 120 piglets on a single filling.
- A complete set includes a docking instrument and a can of lighter fluid in a sturdy cardboard box.

Temperature: approx. 400 °C

Operating time: 1-1.5 hours at maximum temperature

Refilling

It is best to always fill the cylinder before use to eliminate the chance of running out of gas midway through the work. Under normal conditions of use, at least 10 litters of piglets can be docked with one filling. Refilling is very easy: just hold the can of gas upside down on the docking instrument and depress the valve until the gas begins to flow.

The docking instrument is designed to dock the tails quickly and painlessly.

Light the docking instrument with a match or lighter before you enter the birthing pen. Before lighting, the red button is pushed forward to the 'open' position and the gas turned to 'HI' (high). Hold the match or lighter under the wire mesh on the knife until the gas starts to burn.

First the gas burns outside the wire mesh, but after 10-30 seconds the flame catches inside, where it can be observed as an orange glow. While the flame burns from the outside to the inside of the wire mesh, it is best to hold the docking instrument in your right hand, resting on a table.

It will now take 1 to 2 minutes before the knife reaches the right temperature (400 °C) and is ready for use. To save gas, the docking instrument can be turned to 'LO' (low) between each litter. When you are ready to start the next litter of piglets, turn it back to 'HI' and the instrument is immediately ready for use.

When docking tails, the heat of the knife does the actual work. Squeeze gently and carefully until the tail falls off. Do not force the knife through the tail, as this may result in bleeding and can lead to infection.

By using the docking instrument you do not run the risk of damaging the nerves in the tail or tearing them off. The pig retains feeling in the tail, which eliminates the problem of tail biting.

When you are finished using the docking instrument, push the red button back to close off the gas supply.

Replacing the blade

1. Loosen the two gas cylinder attachment screws.
2. Remove the gas cylinder.
3. Note: Do not remove the two o-rings that properly position the gas cylinder.
4. Loosen the blade using a cloth and try not to use any tools; otherwise you may damage the instrument.
5. Screw the new blade in place while ensuring that the sharp knife-edge is aligned with the middle of the red open/close knob. Tighten it securely. Replace the two o-rings behind the red open/close knob. Place them so far back that they fall within the metal housing.
6. Place the gas cylinder back in the metal housing, ensuring that the knife blade is directly above the pressure plate and the red open/close knob can be operated easily.

7. Tighten the two screws, so the gas cylinder cannot turn in the housing.

WARNING

The catalyst inside the knife ensures good, uniform heating; it must not be adjusted. The instrument must not be cooled with water. If it is necessary to clean the knife, use a dry cloth or fine emery cloth if necessary.

The gas cylinder contains gas under pressure and therefore must not be exposed to temperatures above 50 °C or excessive sunlight. If the outside flame does not disappear within 30 seconds, you may need to replace the knife, because the catalyst may be damaged.

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© Copyright 2014