MS Pulse Dosing Unit

1406500, 1406501, 1406502, 1406503
1406504, 1406505

NL	MS Pulse Doseerunit 50 & 250	
	Gebruiksaanwijzing	3
EN	MS Pulse Doseerunit 50 & 250	
	Manual.	31
DE	MS Pulse Doseerunit 50 & 250	
	Gebrauchsanweisung	59
FR	MS Pulse Doseerunit 50 & 250	
	Mode d'emploi.	87
ES	MS Pulse Doseerunit 50 & 250	
	Instrucciones para el uso	115

MS Pulse 50, MS Pulse 250

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MS Pulse 50, MS Pulse 250

1.	General safety instructions.....	32
2.	Introduction	33
3.	Product specifications	34
4.	Pulse 50 and Pulse 250 needle-free injection system.....	35
4.1	Recommended injection site	35
4.2	Cleaning procedure	36
4.3	Use of the Pulse 250.....	37
4.4	Use of the Pulse 50.....	38
4.5	Pressure settings	39
4.6	Injection gun	40
4.6.1	Manual activation and depressurising of the injection gun	40
4.6.2	Contingency procedure: depressurising the injection gun manually	41
4.6.3	Injection gun maintenance (P250).....	41
4.6.4	Injection gun maintenance (replacing O-ring on P50).....	45
4.6.5	Parts of the inlet valve (P50)	49
4.6.6	Parts of the non-return valve (P50)	50
4.6.7	Cleaning the nozzle.....	51
4.6.8	Opening the injection gun	52
4.6.9	Closing the injection gun	53
4.7	Dosing unit	56
4.7.1	Dosing unit specifications.....	56
4.7.2	Dosing unit maintenance	56

1. General safety instructions

This manual is intended to support the use and maintenance of the various systems on the MS Multi treatment unit. Please read and understand the entire manual before using any of these systems.

The safe operation of equipment on the treatment unit is the responsibility of the user. All safety instructions must be complied with before using the various devices on the treatment unit.

NEVER allow untrained personnel to operate devices. **NEVER** leave devices unattended.

The Multi treatment unit and its accompanying devices have been developed for use on animals; **NEVER** use for other purposes.

Take appropriate precautions (e.g. protective clothing, safety goggles, dust mask, ear defenders, gloves).

Pulse 50 and Pulse 250:

NEVER inject liquid into the air. **ALWAYS** point the injection gun in a safe direction. **NEVER** aim at another person.

ALWAYS set the On/Off air supply valve on the dosing unit to OFF when the MS Pulse injection system is not in use. The high-pressure hose and injection gun should also be depressurised.

NEVER store at temperatures **below 0 °C**.

TAKE CARE The stream of liquid from the P50 and P250 injection gun is powerful enough to penetrate gloves, boots and clothing.

NB: The Schipmatic Pulse 250 has been developed to administer mass injections of iron preparations, commonly used vaccines and serums. ***Certain antibiotics and high-viscosity liquids may cause technical problems. Injection by MS Pulse 250 is not recommended in these cases.***

2. Introduction

Congratulations! You are the proud owner of an MS Pulse. The revolutionary systems on the treatment unit were specially developed for use in the pig sector.

Advantages of working with the Multi treatment unit:

1. Hygienic: eliminates cross-contamination by needles or treatment crates,
2. Efficient: significantly reduces the time taken to process a litter of piglets,
3. Animal friendly: reduces pain and stress for the animals.

The results are better growth, fewer losses and lower antibiotic use.

The Schipmatic Multi treatment unit can be equipped with the following modules:

- TagMatic ear tagger
- Needle-free Pulse 250 injection system
- Smart tail docker
- Needle-free P50 injection system
- Cable reel
- Hand hygiene module

This manual explains the use and maintenance of the various modules.

3. Product specifications

Needle-free Pulse 50 and Pulse 250 injection system:

DOSING:	Dosing options for the P50: 0.1-0.5 ml (in increments of 0.0125 ml). Dosing options for the P250: 0.5, 0.75, 1.0, 1.25, 1.50, 1.75, 2.0, 2.50 ml.
INJECTION:	intramuscular/intradermal.
INJECTION TIME:	P50: under 200 milliseconds for 0.1 to 0.5 ml injection. P250: under 200 milliseconds for 0.5 to 2.5 ml injection.
NOZZLE:	Diameter standard 0.20 mm for P250 iron injection, P50 and P250 vaccines standard 0.27 mm.
POWER SOURCE:	The MS Pulse needle-free system operates entirely on compressed air.

4. Pulse 50 and Pulse 250 needle-free injection system

The use of these two systems is virtually identical. The following chapter explains the use and maintenance of Pulse needle-free systems.

In this text, Pulse 50 and Pulse 250 are often abbreviated to P50 and P250.

4.1 Recommended injection site

The recommended injection site is in the neck behind the ear, ideally inside the triangle indicated (not too close to the ear).

Make sure that the safety catch is OFF and that the trigger lock of the injection gun is kept constantly depressed during injection. The injection is activated by pressing the gun against the pig with sufficient force. **Make sure that you always place the injection gun flat against the pig!**

If there are any scars, dirt, abscesses, etc. at the injection site shown, we recommend injecting on the other side of the pig. If there are injuries or damage here too, inject at least 3 cm away from these.

Depending on the liquid used, a swelling or rash may develop at the injection site.

Press the pig against the injection gun with sufficient force to activate the gun. Do not withdraw the pig too quickly, or the injection process will not be completed.

Make sure the correct pressure is set on the dosing unit. If the pressure is too low, the injection will not be completed or liquid may leak from the gun after injection; excess pressure may result in tissue damage or liquid return.

NEVER drop the nozzle or injection gun, or knock them against a hard surface; damage may cause equipment to malfunction. Do not attempt to remove blockages using a needle or other hard object.

Make sure there is no dirt in the inlet valve as this can lead to blockages.

Check the injection visually. The injection pressure may need to be adjusted (depending on the weight and condition of the pig).

4.2 Cleaning procedure

Hygiene in and around needle-free treatment equipment

Pay close attention to hygiene in and around the equipment. Needle-free equipment should be treated in the same way as traditional syringes, so thorough cleaning after use is important.

Follow these steps closely to guarantee even better hygiene of your needle-free system.

Protocol for flushing and cleaning needle-free treatment equipment


1. Start by cleaning the outside of the equipment thoroughly with hot water and detergent (do not submerge).
2. Attach a bottle of MS Rinse Liquid to the pump unit. Pump until the liquid emerging from the nozzle is completely clean.
3. Then attach a bottle of MS Contact Liquid to the pump unit. Pump until liquid emerges from the nozzle (the time taken depends on the length of the dosing hose; see diagram below).
4. Leave this liquid in the nozzle until the next treatment. Do not remove the bottle of MS Contact Liquid from the pump unit (tests have shown that MS Contact Liquid does not have an adverse effect on the nozzle parts).
5. Disinfect all outside parts using an alcohol-based disinfectant.
6. Before starting the next treatment, attach a bottle of MS Rinse Liquid to the pump unit and pump until the MS Contact Liquid has been flushed from the system. To be sure that all liquid has been flushed from the system, use the figures in the diagram below. Take account of the dosing plug fitted.
7. Attach the product to be injected to the pump unit and pump until it emerges from the nozzle. To be sure that all MS Rinse Liquid has been flushed from the system, use the figures in the diagram below. Check the pressure setting, and the system is ready for use.

STORAGE RECOMMENDATION: if the Schipmatic Pulse is being stored away and will not be used for more than ten days, we recommend removing the injection pin from the injection gun and the dosing pin from the dosing units (this extends O-ring lifespan). See instructions above.

NEVER store at temperatures below 0 °C.


4.3 Use of the Pulse 250

1. Attach the high-pressure hose to the output valve as shown opposite. Spray the coupling with silicone oil first (makes assembly easier and extends O-ring lifespan). Tighten firmly (do not use tools). Place the other end of the high-pressure hose in an empty bottle or similar.

2. Attach the bottle as shown above (do not use tools).
3. Attach the air hose to the compressor and the dosing unit.

4. Set the On/Off valve to the On position, holding the high-pressure hose in the empty bottle. The dosing unit will now start to repeat automatically. Liquid will emerge from the hose after a period of time depending on the dose setting and the length of the hose. The liquid will be a mixture of vaccine and water at first but very soon only vaccine will emerge. When this happens, set the valve to Off.

5. Connect the other end of the high-pressure hose to the injection gun. Tighten by hand (do not use tools). The hose should also be disconnected by hand. If it cannot be disconnected, there may still be pressure in the hose. Depressurise by closing the valve and activating the injection gun manually (see instructions above).

6. Set the On/Off valve to the On position. Always aim the injection gun in a safe direction.
7. Adjust the pressure on the dosing unit if necessary. Pull out the yellow knob and turn until the desired pressure has been set. Then push the yellow knob back in.

8. Activate the injection gun manually.
9. Check that the pressure is still at the desired setting.

4.4 Use of the Pulse 50

Setting the dose:

The P50 has a 'dosing knob' on the front of the gun. This can be used to set the dose between 0.1 cc and 0.5 cc.

Turning the dosing knob moves the gun's dosing chamber in or out. The corresponding dose setting is marked on the dosing chamber.

In the picture opposite, the P50 is set to 0.3 cc.


Preparing the P50 for use:

Activate the P50 manually two or three times; it will then draw up the liquid automatically. The P50 is ready for use. (See section 4.6.1 for how to activate the P50 manually.)

If the P50 does not draw up liquid automatically, remove the bottle assembly part from the gun but leave the inlet valve in place. Do not lose the steel ball-bearing in the inlet valve!

Fill the inlet valve with water.

Fit the bottle assembly back on the gun. Activate the gun manually twice.

4.5 Pressure settings

piglet injection gun
(blue or hands-free)

P250

Type of pig	Pressure (bar)
Newborn	
5-15 kg	1.7-2.2
15-25 kg	2.5-3.0
	3.0-3.5


Type of pig	Pressure (bar)
Newborn	
5-15 kg	4.5-5.0
15-25 kg	5.0-5.5
	5.5-6.0

Note: settings for high-viscosity liquids can be approx. 0.5 bar higher.

4.6 Injection gun


The following paragraphs explain the use and maintenance of the P250 injection gun. All of this information also applies to the hands-free variant. Its internal mechanism is exactly the same as in the injection gun.

Similar actions for the Pulse 50 are also explained.

NOTE: viewed from left side of injection gun (side with the RED ring).

4.6.1 Manual activation and depressurising of the injection gun

Manual activation of the injection gun requires two hands. **Be sure to point the handpiece in a safe direction.** Make sure the safety catch is in the off position. Hold the trigger lock in using your index finger. With your other index finger, draw the supply pipe back (see yellow arrow in picture above). The injection gun will now be activated. **Watch out for the stream of liquid.** Take your finger off the supply pipe to allow it to return to its original position.

Once the dosing unit has reloaded, the injection gun is ready for the next manual activation.

Warning: Keep fingers clear of the nozzle!

Manual activation of the Pulse 50 is similar to that of the Pulse 250.

Point the P50 in a safe direction, preferably into a collecting pan.

To activate the P50 once, use your index finger to draw the mechanism back (see picture).

Watch out for the stream of liquid and keep fingers clear of the nozzle.

Then allow the mechanism to spring back to its original position.

The P50 is now ready for the next activation.

4.6.2 Contingency procedure: depressurising the injection gun manually

If the **nozzle is blocked**, for example, you may not be able to depressurise the injection gun in the way described above. If the nozzle is undamaged, it can be disassembled and cleaned. To do so, the injection gun must first be depressurised. (The system must also be depressurised in order to disconnect the high-pressure hose.) To depressurise, attach the supplied spanners to the injection gun and nozzle and give the nozzle approximately one turn to unscrew it. This will cause liquid to drip out slowly. Wait until this has stopped and the system is depressurised. Wear protective clothing and face protection.


4.6.3 Injection gun maintenance (P250)


Replace O-rings after every 4,000 injections or once a month.

Make sure the injection gun is clean and is not connected. Place the handpiece on a clean, flat and stable surface to minimise the risk of losing parts. Make sure that parts cannot become damaged, as this can lead to malfunctioning of the injection gun.

Parts of the injection gun

1. Use the two supplied spanners to open the injection gun.

2. Remove the end piece.

3. Remove the injection pin housing from the injection gun using the 'fork' tool supplied. **Do not use pliers or similar.** Make sure not to damage the housing or injection pin.

If the injection pin housing is difficult to remove, there is probably dirt in the injection gun. Remove any dirt.

4. Push the injection pin out of the housing using the tool supplied.

5. Remove the O-ring on the injection pin using the tool supplied. **Take care: do not damage or remove the white spacer rings.**

6. Fit the new O-ring as shown below, using the tool supplied.

Make sure the O-ring is in the correct position!

Needle housing

7. Replace the O-rings on the needle housing as shown below.

8. Once the O-rings have been replaced, wet the injection pin with water or Medispray and then insert into the housing as shown below.

9. Replace the O-ring on the end piece if it is visibly damaged. Wet the end piece and insert it into the injection pin as shown below.

10. Then insert the assembled needle housing into the injection gun. Wet the assembly with water before insertion.

11. Re-attach the nozzle and nozzle nut as shown below. Make sure that the O-ring on the end piece is present and that the nozzle sits in the nozzle nut correctly.


12. Tighten the nozzle nut using the two spanners (do not use the dosing pipe).


Ready for use. To remove any dirt or contamination, inject with sterile water several times before using vaccines.

4.6.4 Injection gun maintenance (replacing O-ring on P50)

Make sure the gun has been cleaned according to the cleaning procedure and is not connected. When carrying out maintenance on the gun, always work on a clean, flat and stable surface to minimise the risk of losing parts.


1. To remove the dosing chamber from the gun, turn the dose-setting knob clockwise. This extends the dosing chamber to the maximum length.

Once the dosage chamber has been disassembled, the piston can be removed. This is done by sliding the piston out from the 'horseshoe'.

Base of the piston


'Horseshoe' inside the P50

2. The O-ring can now be removed using pincers or another tool.


Take care: do not damage or remove the white spacer rings.


3. The new O-ring can now be slid into place over the piston.

4. Make sure the white supporting ring sits behind the O-ring as shown below.

5. The P50 can now be reassembled. The first step is to put the piston back into the P50. The piston must be fitted back into the horseshoe. It is important that the piston should sit securely in the horseshoe as shown below. If this is not the case, serious damage may be caused to the P50.

6. The dosing chamber can now be returned to the P50. Make sure the O-ring has been wetted with water!


7. Turn the dose-setting knob so that the red dot is at the top.

8. Slide the dosing chamber over the piston and give it a half turn clockwise. The dosing chamber is now back in upright position as shown below.


9. The dosing knob can now be used to set the right dose and the Pulse 50 is ready for use again.

4.6.5 Parts of the inlet valve (P50)

NB: The inlet valve has grooved parts to prevent the bottle assembly from turning during use of the P50. Make sure that these parts fit together correctly when the bottle assembly nut is tightened, or the inlet valve will not seal properly.


Make sure that all parts are coated thoroughly with silicone spray during assembly. This extends the lifespan of the O-rings and helps the parts to slide together more easily.

TAKE CARE not to lose the steel ball-bearing. A spare ball-bearing is included with the set of O-rings for the P50.

4.6.6 Parts of the non-return valve (P50)

Parts of the non-return valve:

The non-return valve requires little or no maintenance. However, if a blockage does occur, try to clear it first by cleaning the nozzle (see section 4.6.7).

Only when all obvious causes of blockage have been ruled out should you disassemble and clean the non-return valve. TAKE CARE not to lose any of the small parts.

Coat all parts thoroughly with silicone spray before assembly. Use the supplied spanners to tighten the nozzle nut.

4.6.7 Cleaning the nozzle

Dirt normally collects in the inlet valve filter. This stops the nozzle from becoming blocked.

Nozzle P50 + P250

Inlet valve filter P250

The different injection guns have nozzles with different aperture diameters. If a nozzle does become blocked, follow the cleaning procedure described below.

Insert the nozzle into the injection gun the other way round and activate the gun manually several times. Increase the pressure if necessary!

Nozzle in correct position


Nozzle back-to-front for cleaning

4.6.8 Opening the injection gun

It may be necessary to open the handpiece because the metal core is no longer moving freely within the plastic shell. This may occur due to use in extremely dirty conditions, poor cleaning, no cleaning, or use of "sticky" products such as iron.

1. Place the handpiece on a clean, flat and stable surface to minimise the risk of losing parts.
2. Slide the rubber cover (and O-ring) towards the nozzle.
3. Remove the five screws (note that they are not all the same length).

4. Turn the handpiece over and lay it on the surface (with the side that held the screws facing down). Carefully separate the two halves. Watch out: the spring, bush and trigger may spring out of the shell. Take care not to lose these.


5. The core can now be removed from the shell. Vaccine residues on the core can now be removed with soap, water and paper towels. **DO NOT submerge the core in water.**


6. Clean the shells with soap, water and paper towels. It is OK to submerge shells. Take care not to lose any of the small parts.
7. Dry the parts thoroughly. The plastic shells need not be completely dry.

4.6.9 Closing the injection gun




1. Place the right-hand shell on the surface. This is the side with the holes for the screw heads.
2. Insert the spring and bush into the plastic shell.

Position the trigger lock as shown.


3. Hold the lock in this position with one hand; with your other hand, move the spring to the position shown below. This can be difficult and may take a number of attempts.


The correct spring position is as shown in the picture.

4. Now fit the safety catch as shown below. Make sure that the trigger lock and spring remain in position.


5. Position the core in the shell. Make sure that the trigger lock and spring remain in position.


6. Fit the other half of the shell. Make sure that the trigger lock and spring remain in position.

7. Test the operation of the trigger lock.

8. If the lock cannot be pushed in completely, the spring may not be correctly positioned. If this is the case, open the handpiece and start again from step 1.
9. Turn the handpiece over and insert the five screws (note the different lengths). Do not over-tighten the screws.

10. Slide the rubber cover and O-ring back into position (see below).

Groove


O-ring in groove

4.7 Dosing unit

4.7.1 Dosing unit specifications


P250 Dosing chamber

4.7.2 Dosing unit maintenance


Replace O-ring after every 4,000 injections or once a month

1. Make sure the dosing unit is clean and is not connected. Place it on a clean, flat and stable surface to minimise the risk of losing parts. Make sure that parts cannot become damaged, as this can lead to malfunctioning of the injection gun.
2. Open the dosing chamber as shown below, using the spanner supplied.


3. Pull the dosing chamber around 10 mm upwards. Remove the dosing pin from the piston (horseshoe shape). See picture above. Now remove the whole dosing chamber from the dosing unit.


- Remove the piston from the dosing chamber.

- Remove the O-ring on the dosing pin. **Take care: do not remove or damage the white spacer rings.**
- Once the O-ring has been removed, wet the dosing pin with water and insert it into the dosing chamber. Note that the second ring should still be just visible.


- Insert the dosing chamber into the dosing unit as shown below.

- The dosing pin should be fitted into the horseshoe as shown below. Pull gently on the dosing chamber; you should feel resistance.

9. When you are sure that the dosing chamber is sitting securely in the horseshoe, push the dosing chamber downwards. Now pull the dosing chamber upwards again, by around 10 mm, to check that it is sitting securely in the horseshoe. Now push the dosing chamber tightly into the dosing unit.

10. Tighten the nut (hand-tight) using the spanner supplied.

Note: to prevent damage or malfunctioning of the dosing unit, it is important to fit the dosing pin securely into the horseshoe.

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**Schippers Italia Srl**
**Schippers Schweiz GmbH**
**Schippers UK Ltd.**
**LLC Schippers Russia**
**Schippers Canada Ltd**
7102 52nd Street Bay # 18 • Lacombe, AB • T4L 1Y9