Reason for use:

Ketosis occurs in dairy cows as a consequence of the negative energy balance in the early postpartum period, and can result in lowered milk production, impaired fertility, and increased risk of periparturient diseases, such as displaced abomasum. Subclinical ketosis usually precedes clinical ketosis, and is much more common. Subclinical ketosis causes greater economic losses on a herd basis than clinical ketosis. Subclinical ketosis can be diagnosed by measuring ketone bodies present in blood, urine, or milk. The Keto-Test™ (ketolac test strip) is a convenient, very effective method of detecting cows with subclinical ketosis.^{1,2}

Test principle

Beta-hydroxybutyric acid in the milk passes through the reagent part of the test strip and is converted by beta-hydroxybutyric acid dehydrogenase (BHBHD) to acetoacetic acid (AcAc). The NADH produced from NAD in the process reduces nitroretazolum blue (NTB) to formazan, which is purple in colour. The concentration of beta-hydroxybutyric acid in the milk can then be estimated from the change in colour.

Directions for use:

1. Remove test strip container from the refrigerator and bring to room temperature. Immediately close the container again firmly.
2. Dip a single test strip for 3 seconds into a small quantity of fresh milk (in a clean container). If the milk sample has been standing for some time it must be thoroughly shaken before carrying out the test. It, too, should be at room temperature.
3. Shake twice vigorously to remove excess milk.
4. After 1 minute read test strip color and compare to colour chart on bottle label (if colour reaction varies at the margin, base the reading on the central colour).

Interpretation:

Milk BHBA concentration

0-50 µmol/l
100 µmol/l
200 µmol/l
>500 µmol/l

Keto-Test™ Assessment

Normal (-)
Weak positive (+/-)
Positive (+)
Highly positive (++)

This milk test strip system semi quantitatively measures Beta-hydroxybutyric acid (BHBA), which is a key ketone body present in milk. The concentration of BHBA present in the milk can be estimated by the strength of the colour change reaction. The more purple the result, the higher the concentration of ketones present in the milk and the greater the likelihood of ketosis. BHBA levels > 200 µmol/l of milk were 4 times more likely to come from cows with subclinical ketosis. This Keto-Test™ was found to be the best cowside method for measuring BHBA.^{1,2}

Cautions:

- This test is designed for measuring ketones present in milk. If other fluids are measured, then test result interpretations may be different.
- If milk cannot be tested immediately following collection, store in a refrigerator and bring to room temperature before testing.
- If the sample being tested is colostrum, the reading of the colour reaction is more difficult because of the milk's intrinsic colour.

Storage and handling:

- Store at 2 - 8°C (36 - 46°F).
- The product shelf life is 1.5 years from date of manufacture.
- Continuous exposure of the product to temperatures above 25°C (77°F) results in a shelf life of only 84 days.
- Keep tubes of strips cool until use.
- Do not use test strips that have discoloured after extended storage.
- Keep the test strip container tightly closed.
- Do not expose test strips to direct sunlight.
- Do not touch the test strip reagent area.

1. Geishauser T, Leslie K, Tenhag J and Bashiri. Evaluation of Eight Cow-side Ketone Tests in Milk for Detection of Subclinical Ketosis in Dairy Cows. 2000. *J. Dairy Sci.* 83:298-299.
2. Geishauser T, Leslie K, Kelton D and Duffield T. Evaluation of Five Cow-side Tests for Use with Milk to Detect Subclinical Ketosis in Dairy Cows. 1998. *J. Dairy Sci.* 81:438-443.

Manufacturer:
SANWA KAGAKU
KENKYUSHO CO. LTD
Nagoya, Japan

For additional product information call + 44 (0)1256 353131 or email: elancouk@lilly.com

Elanco, Keto-Test and the diagonal bar are trademarks owned or licensed by Eli Lilly and Company, its subsidiaries or affiliates.

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