MastDecide – Decision-making support for evidence-based mastitis treatment

Objective

This quick test makes it possible to check the sterility of quarter milk samples taken from dairy cows with clinical mastitis. This test is not intended to replace cytomicrobiological examinations in a laboratory.

Contents

Ready-to-use test with liquid test culture media (agar) in containers (medium 1: white cap; medium 2: yellow cap), single-use pipettes, autoclavable waste bag, sterile milk sampling test tubes (for quarter milk samples), and udder/teat cleaning wipes.

Other required materials (not provided)

Disposable paper towels and disinfectant (70% ethanol) for teat cleaning and disinfection, disinfectant for surface cleaning of working area, disposable gloves, sample rack, incubator (37 °C $\pm$ 2 °C).

Warning

The product contains irritant substances in highly diluted concentrations. Avoid direct contact with the skin. Wear gloves when using the product. In case of accidental contact with the eyes, rinse the affected eye gently with tap water for several minutes, remove contact lenses if worn, and rinse again. Consult an optician/ophthalmologist if necessary.

Used test-media containers must not be reopened and should be disposed of appropriately (rendered harmless prior to disposal with regular waste) in order to prevent environmental contamination with the micro-organisms that they contain.

Not for human consumption. Keep out of the reach of children.

Instruction manual

Milk sampling

Milk samples should be taken in as sterile a manner as possible in order to avoid false or inaccurate test results.

We recommend the following procedure:

- Wear disposable gloves during sampling.
- Clean the udder and teats in as dry a manner as possible to remove visible dirt.
- Discard the first three jets of milk per teat.
- Next, disinfect the teat, especially around the teat apex and the teat canal opening, using a disposable cloth soaked in a suitable disinfectant.
- Allow the teats to dry prior to sampling.
- During sampling, hold the milk sampling container as horizontally as possible. Opening and closing of the container should take place under the teat. Avoid touching the inside of the cap of the milk sample container (either with your hands or the teat).

Testing

Attention

Wear clean disposable gloves during testing. A clean working environment and the effective disinfection of the work surface before and after using the test materials is also recommended. The working area should be separated from rooms where food for humans is stored or consumed.

Preparation

Switch on the incubator. The working temperature is 37 °C $\pm$ 2 °C.

1. Use two test containers (test medium 1: white cap, test medium 2: yellow cap) for each sample (see steps 3. to 7.).

2. Gently mix (swish) the quarter milk sample that was obtained under aseptic conditions.

3. Under aseptic conditions and using the single-use pipette, take 0.1 ml of the well-mixed milk sample. The pipette must be completely full of milk.

4.

4. Carefully open the first test medium container (test medium 1: white cap). Do not touch the inner surfaces, only the outside of the container and lid.

5. Add the 0.1 ml milk sample to test medium 1 (white cap). Avoid touching the inside of the container with the tip of the pipette.

6. Carefully close the cap of the test medium and gently mix it (by swishing). Do not put down the pipette due to the risk of contamination. The test medium must not be reopened. Incubate at 37 °C.

7. Repeat steps 3. to 6. with test medium 2 (yellow cap).

8. Label the caps as required.

Analysis of the test samples


The test samples are examined after an incubation period of 12 to 14 hours. Assess the colour of both test media. Use the table below to assess the sterility of the samples tested:

Test medium 1 – white cap (general pathogen growth)	Test medium 2 – yellow cap (Gram-positive bacteria)	Test result
white	pink	Gram-negative/Coliform ^a
white	white	Gram-positive ^b
pink	pink	No pathogens detected

^a *E. coli*, *Klebsiella* spp., other coliform counts (*Citrobacter*, *Enterobacter*)

^b *S. aureus*, CNS, *S. agalactiae*, *S. dysgalactiae*, *S. uberis*, enterococci, *Lactococcus lactis*

No growth of: *Corynebacteria*, *Pseudomonads*, *Trueperella pyogenes*, yeasts, mycoplasmas

Notes of assessment

The test medium only changes colour from pink to white if microorganisms are growing in it. The change in colour might not affect the whole of the test medium; even if only parts have changed colour, this also indicates micro-organism growth and the result is considered to be positive. If both test media stay pink, no pathogen growth is detectable. If both test media change colour from pink to white, gram-positive microorganisms are present. If only test medium 1 (with the white cap) changes colour, gram-negative pathogens are likely to be present. In this case, we recommend incubating the test samples for a further two hours before re-evaluating (to increase the test sensitivity). If test medium 2 stays pink, gram-negative pathogens are present; if it turns white, gram-positive pathogens are present.

Storage and shelf life

The fresh test media are transparent with a distinct pink colour. If the fresh, unused media have changed colour or have become cloudy, they should not be used and should be disposed of.

The test media should be stored in a dark and cool place (preferably between 4 and 8 °C). Don't freeze them. The shelf life when stored correctly is up to 9 months. Do not use the test materials after the expiry date. The test media are intended for single use only. Keep used test media out of the reach of children.

Disposal

After assessing the test samples, the test media must be rendered harmless and disposed of safely. For intermediate storage, please place them in an autoclavable disposal bag. We recommend heating to 121 °C with at least 15 minutes at that temperature in a pressured atmosphere (e.g. in an autoclave). The test media that have been rendered harmless may be disposed of with regular household waste. Alternatively, the used test media may be disposed of by a veterinarian.

Troubleshooting

1. The test shows negative results or gram-negative bacterial growth, in contrast to gram-positive bacterial growth results from an official laboratory:

Please check the following points:

- Please pay attention to ensuring a high level of hygiene during sampling. Any contamination may lead to a false positive result, as bacteria might interfere with each other and stop each other from growing.
- The incubator should be pre-set and warmed up to 37 °C when placing the test agars.
- Check the incubation temperature with a thermometer, ensuring that it is at 37°C. If not, please correct.
- When adding milk to the test agar, check whether air was drawn into the pipette. If so, please repeat procedure, as there may be an insufficient amount of test media in the agar, leading to a negative test result.
- Only evaluate the test after 12 hours: an earlier reading may lead to an increased number of gram-negative and negative results.
- If an agar turns white at certain zones, it is regarded as positive (please note that not all the agar has to turn white for a result to be regarded positive).
- Don't shake or mix the test agar immediately prior to evaluation, as this may lead to a higher likelihood of more gram-negative growth or no growth.
- If the results show 'Gram-negative', you may incubate for another 2 hours (total of 14 hours) and re-evaluate the results. If so, record the results under 'notes'

2. The test shows gram-positive or gram-negative results, but no bacteria have been detected in a laboratory. Below are the possible reasons that explain this discrepancy:

- In laboratory assays, only 10 µl of milk are usually evaluated, whereas 100 µl are evaluated in this quick test. In light of this, the quick test may show up as correctly gram-positive or gram-negative in some of the samples that were tested negative in the laboratory.
- Possibly contaminated samples may be evaluated in the laboratory as having no growth or relevance, whereas in the quick test this contamination may lead to gram-positive or gram-negative results.

Further information

All information above was provided in good faith and to the best of our knowledge. None of the information contained herein may be interpreted as a guarantee or warranty. 'As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the veterinarian after all clinical and laboratory findings have been evaluated'.

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Veterinary Enterprises Europe B.V.